



35543

## Aluminium Hydroxide Gel extrapure AR, ExiPlus, Multi-Compendial

Part A

CAS : 21645-51-2

Molecular Formula : Al (OH)<sub>3</sub>

Molecular Weight : 78.0

### Specifications

ExiPlus	meets compendial specs of IP
Appearance (Colour)	White to almost white
Appearance (Form)	Light amorphous powder
Assay	47 - 60% (as Al <sub>2</sub> O <sub>3</sub> )
pH (4% aq. suspension)	max.10.0

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	28183000 (GST 18%)
500 Gms	28183000 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

5 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com