



35889

Potassium Tellurite Monohydrate extrapure, 90%

Part A

CAS : 123333-66-4

Molecular Formula : K₂O₃Te.H₂O

Molecular Weight : 271.8 (253.80 anhy)

Specifications

| | |
|---------------------|------------------------|
| Appearance (Colour) | White to grayish white |
| Appearance (Form) | Crystalline powder |
| Assay | min. 90% |

General Information

| | |
|---------------------|---|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:3284 , IMCO Class No.:6.1 , Packing Group:III |
| HSN Code | |
| 10 Gms | 28429010 (GST 18%) |
| 100 Gms | 28429010 (GST 18%) |
| 25 Gms | 28429010 (GST 18%) |
| Type of Packing | |
| 100 Gms | Glass Bottle |
| 25 Gms | Glass Bottle |

Available Packages

10 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com