



35981

L-Prolinamide extrapure, 98%

Part B

CAS : 7531-52-4

Molecular Formula : C₅H₁₀N₂O

Molecular Weight : 114.15

Specifications

Appearance (Colour)	White to yellow
Appearance (Form)	Powder
Assay (NT)	min. 98.0%
Water (KF)	max. 1.0%
Melting Point	98 - 102°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
100 Gms	29333990 (GST 18%)
25 Gms	29333990 (GST 18%)
5 Gms	29333990 (GST 18%)
500 Gms	29333990 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

5 Gms

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com