



36078

## DWCNT Type 1 Carbon Nanotubes Double Walled, Length 50µm, OD 2-4nm, 60%

Part C

CAS : 308068-56-6

### Specifications

Appearance (Colour)	Black
Appearance (Form)	Powder
Assay	min. 60%
OD	2- 4 nm
Length	~ 50µm
SSA	>350m <sup>2</sup> /g

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1325 , IMCO Class No.:4.1 , Packing Group:III
HSN Code	
100 Mg	28030090 (GST 18%)
500 Mg	28030090 (GST 18%)

### Available Packages

100 Mg

500 Mg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com