



36179

## Sorbitan Monostearate (SPAN 60®) extrapure

Part A

CAS : 1338-41-6

Molecular Formula : C<sub>24</sub>H<sub>46</sub>O<sub>6</sub>

Molecular Weight : 430.63

### Specifications

Appearance (Colour)	Pale yellow
Appearance (Form)	Flakes
Melting Point	54 - 57°C
Acid value	max. 10
Saponification value	147 - 157
Hydroxyl value	235 - 260

### Other Information

Notes

Span is a registered trademark of Croda International PLC

SPAN 60 <sup>™</sup> is used as a non-ionic detergent.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	34024200 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com