



36324

Glycerol-3-Phosphate Dehydrogenase ex. Rabbit Muscle, 15U/mg

Part B

CAS : 9075-65-4

Specifications

Appearance (Colour)	White
Appearance (Form)	Hygroscopic crystalline Powder
Activity	min. 15 units/mg powder
Unit Definition	One unit will reduce 1.0 micromole of 3-phosphoglycerate to D-glyceraldehyde-3-phosphate per minute in a coupled system with 3-phosphoglyceric phosphokinase at pH 7.6 at 25°C.
Enzyme solution for activity determination	Prepared in 50 mM BICINE buffer containing 0.1% BSA, pH 9.0.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Units	35079099 (GST 18%)

Available Packages

500 Units

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com