



37163

## Devardas Alloy (45% Al, 50% Cu, 5% Zn) extrapure AR, ACS

Part A

CAS : 8049-11-4

### Specifications

Appearance (Colour)	Grey
Appearance (Form)	Powder
Copper (Cu) content	49-51%
Aluminium (Al) content	44-46%
Zinc (Zn) content	4-6%
Nitrogen (N)	0.005%
Sieve test	min. 95% passing through 200 mesh

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3089 , IMCO Class No.:4.1 , Packing Group:II
HSN Code	
100 Gms	74032900 (GST 18%)
250 Gms	74032900 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
250 Gms	Plastic Bottle

### Available Packages

100 Gms  
250 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com