



37256

Tetrabutylammonium Dihydrogen Phosphate (TBADP) extrapure AR, 98%

Part A

CAS : 5574-97-0

Molecular Formula : C₁₆H₃₈NO₄P

Molecular Weight : 339.45

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline Powder
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless
Assay (T) (on dry basis)	min. 98%
Absorbance (A) of 0.005 M aq. solution in a 1cm cell v/s H ₂ O	
@210nm	max. 0.05
@220nm	max. 0.04
@300nm	max. 0.02
@500nm	max. 0.02

Other Information

Application

Suitable for Ion Pairing & HPLC applications.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29239000 (GST 18%)
5 Gms	29239000 (GST 18%)

Available Packages

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com