



37524

Iodine Resublimed AR, ACS, ExiPlus, Multi-Compendial, 99.8%

Part A

CAS : 7553-56-2

Molecular Formula : I

Atomic Weight : 126.90

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	Bluish black
Appearance (Form)	Crystalline scales/balls
Solubility	Passes test
Assay (Iodometric)	min. 99.8%
Non Volatile Matter	max. 0.005%
Chlorine & Bromine (Cl)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3495 , IMCO Class No.:8(6.1) , Packing Group:III
HSN Code	
100 Gms	28012000 (GST 12%)
500 Gms	28012000 (GST 12%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com