



38516

N,N-Methylene Bisacrylamide (bis-Acrylamide)
3x cryst. extrapure AR, 99.5%

Part A

CAS : 110-26-9

Molecular Formula : C₇H₁₀N₂O₂

Molecular Weight : 154.17

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 2.5% aq. solution	Clear
Solubility (Colour) 2.5% aq. solution	Colourless
Assay (ex N)	min. 99.5%
pH (2.5% aq. solution)	min. 5.0
Specific conductance (2% aq. solution)	max. 10 µmho/cm
Acrylic acid (C ₃ H ₄ O ₂)	max. 0.02%
Absorbance (A) of 1% aq. solution in a 1cm cell v/s H ₂ O	
@290nm	max. 0.2
Water (KF)	max. 0.5%
Particle Size	max. 75 µ

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29241900 (GST 18%)
25 Gms	29241900 (GST 18%)
500 Gms	29241900 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
