



38862

N,N-Methylene Bisacrylamide (bis-Acrylamide) pure, 98%

Part A

CAS : 110-26-9

Molecular Formula : C₇H₁₀N₂O₂

Molecular Weight : 154.17

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 2.5% aq. solution	Clear to Opalescent
Solubility (Colour) 2.5% aq. solution	Colourless
Assay	min. 98%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29241900 (GST 18%)
5 Kg	29241900 (GST 18%)
500 Gms	29241900 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

100 Gms
500 Gms
5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com