



39192

## 2'-Deoxyadenosine-5'-Monophosphate Disodium Salt (dAMP-Na<sub>2</sub>) extrapure, 98%

Part B

CAS : 2922-74-9

Molecular Formula : C<sub>10</sub>H<sub>12</sub>N<sub>5</sub>Na<sub>2</sub>O<sub>6</sub>P

Molecular Weight : 375.19

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min.98%
Spectral Ratio (A250nm/A260nm)	~ 0.80
Spectral Ratio (A280nm/A260nm)	~ 0.14
Spectral Ratio (A290nm/A260nm)	~ 0.007

### General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Mg	29349990 (GST 18%)
Type of Packing	
100 Mg	Glass Bottle

### Available Packages

100 Mg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com