



39227

## tert-Butyldimethylsilyl Chloride (TBDMS.Cl) extrapure, 98%

Part B

CAS : 18162-48-6

Molecular Formula : C<sub>6</sub>H<sub>15</sub>ClSi

Molecular Weight : 150.73

### Specifications

Appearance (Colour)	White
Appearance (Form)	Hygroscopic compound
Solubility (Turbidity) 10% solution in methanol	Clear
Solubility (Colour) 10% solution in methanol	Colourless
Assay	min. 98%
Melting Point	86 - 89°C

### Other Information

Application

Versatile protecting reagent for amines, amides, and especially alcohols, among others.

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2925 , IMCO Class No.:4.1(8) , Packing Group:II
HSN Code	
100 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
500 Gms	29319090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com