



39238

General Purpose Kit - I (Pre-Weighed Media)

Part D

Other Information

Applications

Nutrient agar cultivates a wide variety of microorganism and MacConkey agar is used in isolation of gram negative enteric bacteria.

Composition

(63971)Nutrient Agar	28.00 gm/lt.
(76875)MacConkey Agar	50.00 gm/lt.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com