



39607

## Xanthine Oxidase Suspension ex. Butter Milk, 1-4 U/mg protein

Part B

CAS : 9002-17-9

### Specifications

Appearance (Colour)	Brown
Appearance (Form)	Liquid suspension
Activity	1.0 - 4.0 units/mg protein
Unit Definition	One unit will convert 1.0 micromole of xanthine to uric acid per minute at pH 7.5 at 25°C.
Enzyme solution for activity determination	Prepared in 50 mM potassium phosphate buffer, pH 7.5.

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	36 months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Units	35079099 (GST 18%)
5 Units	35079099 (GST 18%)

### Available Packages

5 Units  
25 Units

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com