



39798

## Resorcinol ACS, 99%

Part B

CAS : 108-46-3

Molecular Formula : C<sub>6</sub>H<sub>6</sub>O<sub>2</sub>

Molecular Weight : 110.11

### Specifications

Appearance (Colour)	White to slight pink
Appearance (Form)	Crystals or Flakes
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Almost colourless
Assay (GC)	min. 99%
Insoluble matter	max. 0.005%
Residue after ignition	max. 0.01%
Melting Point	109 - 112°C
Titration acid	max. 0.004 meq/g

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2876 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	29072100 (GST 18%)
500 Gms	29072100 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

### Available Packages

100 Gms

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com