



40057

Ammonium Phosphomolybdate Hydrate extrapure, 61% Mo

Part B

CAS : 54723-94-3

Molecular Formula : $(\text{NH}_4)_3\text{PMo}_12\text{O}_{40} \cdot x\text{H}_2\text{O}$

Molecular Weight : 1876.35 (anhydrous basis)

Specifications

Appearance (Colour)	Faint yellow to dark yellow
Appearance (Form)	Powder or crystals
Molybdenum (Mo)	~ 61 %

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
100 Gms	28429090 (GST 18%)
25 Gms	28429090 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
25 Gms	Plastic Bottle

Available Packages

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com