



40180

Potassium Hexafluorotitanate pure, 98%

Part A

CAS : 16919-27-0

Molecular Formula : F6K2Ti

Molecular Weight : 240.05

Specifications	
Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Assay	min. 98%

General Information	
Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3288 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	28269000 (GST 18%)
5 Kg	28269000 (GST 18%)
500 Gms	28269000 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages	
100 Gms	
500 Gms	
5 Kg	

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com