



40251

Barium Sulphate extrapure AR, ExiPlus, Multi-Compendial, 98%

Part A

CAS : 7727-43-7

Molecular Formula : BaSO₄

Molecular Weight : 233.40

Specifications

ExiPlus	Meets compendial specs of Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min. 98%
Iron (Fe)	max. 0.005%
Heavy Metals (Pb)	max. 0.005%
Loss on ignition	max. 1%
Acid insoluble matter	max. 0.3%
Acidity	max. 0.1mN%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	28332700 (GST 18%)
500 Gms	28332700 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com