



40474

BOC-D-Phenylalanine extrapure, 99%

Part B

CAS : 18942-49-9

Molecular Formula : C₁₄H₁₉NO₄

Molecular Weight : 265.30

Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Assay	min. 99%
Melting Point	78 - 84°C

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)

Available Packages

10 Gms
25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com