



40555

## Phenylhydrazine Hydrochloride pure, 99%

Part A

CAS : 59-88-1

Molecular Formula : C<sub>6</sub>H<sub>8</sub>N<sub>2</sub>.HCl

Molecular Weight : 144.60

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay (T)	min. 99%
Sulphated Ash	max. 0.1%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2811 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
100 Gms	29215990 (GST 18%)
250 Gms	29215990 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle

### Available Packages

100 Gms  
250 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com