



40653

Sspl (10U/ul)

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Colourless
Appearance (Clarity)	Clear
Activity	One unit is defined as the amount of enzyme required to digest 1 µg of lambda DNA in 1 hour at 37 °C in 50 µL of assay buffer. Incubation of single stranded and double stranded radiolabeled oligonucleotides with 10 units of the enzyme for 4 hours at 37 °C.
DNA degradation	

Other Information

Description	Source: Sphaerotilus species 5'-AAT?ATT-3' 3'-TTA?TAA-5'
Includes	250U Sspl - 25 ?l 10X Reaction Buffer – 1.0ml 500U Sspl - 50 ?l 10X Reaction Buffer – 1.0ml

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 Units	38229090 (GST 12%)
500 Units	38229090 (GST 12%)

Available Packages

250 Units
500 Units

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com