



40707

Potassium Oxalate Monohydrate extrapure AR, ACS, 99.5%

Part A

CAS : 6487-48-5

Molecular Formula : $K_2C_2O_4 \cdot H_2O$

Molecular Weight : 184.24

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (Redox Titration)	min. 99.5%
pH (5% aq. solution)	7.0 - 8.5
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.001%
Copper (Cu)	max. 0.0005%
Nitrogen compounds	max. 0.002%
Sodium (Na)	max. 0.02%
Insoluble matter	max. 0.01%
Ammonium (NH ₄)	max. 0.002%
Substances darkened by hot Sulfuric acid	Passes test
Neutrality	Passes test

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Kg	29171190 (GST 18%)
500 Gms	29171190 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
