



40876

## Eriochrome Black T extrapure AR

Part A

CAS : 1787-61-7

Molecular Formula : C<sub>20</sub>H<sub>12</sub>N<sub>3</sub>NaO<sub>7</sub>S

Molecular Weight : 461.39

### Specifications

Appearance (Colour)	Brown to Black
Appearance (Form)	Powder
Loss on drying	max. 7%
Absorption maxima	610 - 625 nm ( at pH 10 )
Sutability as indicator	Passes test

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	32049000 (GST 18%)
25 Gms	32049000 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

### Available Packages

25 Gms

---

100 Gms

---

500 Gms

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com