



40881

## Chlorogenic Acid extrapure, 98%

Part B

CAS : 327-97-9

Molecular Formula : C<sub>16</sub>H<sub>18</sub>O<sub>9</sub>

Molecular Weight : 354.31

### Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Crystalline powder
Assay	min. 98%
Melting Point	205 - 215°C
Loss on drying	max. 5%

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	98020000 (GST 18%)
100 Mg	98020000 (GST 18%)
5 Gms	98020000 (GST 18%)
Type of Packing	
1 Gms	Glass Bottle
100 Mg	Plastic Bottle
5 Gms	Glass Bottle

### Available Packages

100 Mg

1 Gms

5 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com