



41739

Z-Phenylalanine-p-Nitroanilide extrapure

Part B

CAS : 19647-71-3

Molecular Formula : C23H21N3O5

Molecular Weight : 419.4

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Powder
Assay	min 98%

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29379020 (GST 18%)

Available Packages

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com