



42020

Dodecyltrimethylammonium Chloride (DOTAC)
pure, 98%

Part B

CAS : 112-00-5

Molecular Formula : C15H34ClN

Molecular Weight : 263.89

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless
Assay	min. 98 %
Water (KF)	max. 2 %

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
100 Gms	29239000 (GST 18%)
500 Gms	29239000 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com