



42180

## 2,4-Dichlorophenylboronic Acid extrapure, 98%

Part B

CAS : 68716-47-2

Molecular Formula : C<sub>6</sub>H<sub>5</sub>BCl<sub>2</sub>O<sub>2</sub>

Molecular Weight : 190.82

### Specifications

Appearance (Colour)	Off - white
Appearance (Form)	Powder
Assay (NaOH titration)	min. 98%
Melting Point	240 - 245°C

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)

### Available Packages

1 Gms

5 Gms

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com