



42215

Sodium Caprylate extrapure, 99%

Part A

CAS : 1984-06-1

Molecular Formula : C₈H₁₅O₂Na

Molecular Weight : 166.20

Specifications

Appearance (Colour)	White to almost white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (On dry basis)	min. 99%
pH (10% aq. solution)	8.0 - 10.5
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.001%
Water (KF)	max. 3%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29159099 (GST 18%)
250 Gms	29159099 (GST 18%)
500 Gms	29159099 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
250 Gms	Glass Bottle
500 Gms	Plastic Bottle

Available Packages

100 Gms

250 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com