



42418

Dish Evaporating Porcelain, 40ml, OD 60mm, RB, Deep with Spout

Part L

Specifications

| | |
|--------------------------------|---------------------------------------|
| Appearance (Colour) | White to Off white |
| Appearance (Form) | Evaporating dish |
| Material of Construction (MOC) | Porcelain |
| Evaporating Dish Design Type | Deep with Spout, Round bottom |
| Evaporating Dish Top Outer Dia | 60mm \pm 3mm |
| Evaporating Dish Height | 22mm \pm 5mm |
| Evaporating Dish Weight | 28gm \pm 5gm |
| Max Usable Temperature | 1150°C \pm 50°C |
| Flextural Strength @20°C | Min. 80 Mpa (815 Kg/cm ²) |

Other Information

| | |
|-------------|---|
| Application | A high quality evaporating basin for critical laboratory use. |
|-------------|---|

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 25 to 40°C (Room Temperature) |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 6 Nos | 69091100 (GST 18%) |

Available Packages

6 Nos

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com