



42568

Salicylic Acid extrapure AR, 99.9%

Part A

CAS : 69-72-7

Molecular Formula : C₇H₆O₃

Molecular Weight : 138.12

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in methanol	Clear
Solubility (Colour) 10% solution in methanol	Colourless
Identification (IR Spectrum)	Passes test
Assay (by Neutralization)	min. 99.9%
Melting Point	159° - 161°C
Sulphated Ash	max. 0.02%
Chloride (Cl)	max.0.001%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.0001%
Heavy Metals (Pb)	max. 0.0001%
Copper (Cu)	max. 0.0001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2.5 Kg	29182110 (GST 18%)
500 Gms	29182110 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com