



42577

Copper Metal Turnings, 99.5%, 40-100 mesh

Part A

CAS : 7440-50-8

Atomic Formula : Cu

Atomic Weight : 63.55

Specifications

Appearance (Colour)	Reddish brown
Appearance (Form)	Metallic turning
Solubility (Turbidity) 10% solution in hot con. HNO ₃	Clear
Solubility (Colour) 10% solution in hot con. HNO ₃	Blue
Assay (Iodometric)	min.99.5%
Mesh size	40 - 100 mesh

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3089 , IMCO Class No.:4.1 , Packing Group:II
HSN Code	
100 Gms	74061000 (GST 18%)
500 Gms	74061000 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com