



42780

Tin Metal Powder extrapure AR, 99.9%, ~325 mesh

Part A

CAS : 7440-31-5

Molecular Formula : Sn

Atomic Weight : 118.69

Specifications

Appearance (Colour)	Grey
Appearance (Form)	Fine powder
Assay	min 99%
Iron (Fe)	max 0.02%
Mesh size	80% over 325 mesh (45 micron)

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	80011090 (GST 18%)
250 Gms	80011090 (GST 18%)
Type of Packing	
250 Gms	Glass Bottle

Available Packages

100 Gms
250 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com