



42943

Trichloroethylene Dried, 99.5%, 0.005% water

Part A

CAS : 79-01-6

Molecular Formula : C₂HCl₃

Molecular Weight : 131.39

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 10
Assay	min. 99.5%
Density (g/ml) @ 20°C	1.462-1.464
Refractive Index (20°C)	1.477-1.478
Boiling Range	86-87°C
Non Volatile Matter	max. 0.001%
Acidity (HCl)	max. 0.0005%
Water	max. 0.005%

Other Information

Notes

Packed with Dr. Dried Seal which are designed with corrosion free Crown Cap & with specially developed coloured septa. This enable repetitive needle puncturing easy and securing the consistency of the product from outside atmosphere.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1710 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
500 ml	29032200 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
