



43068

Triphenylphosphinegold(I) Chloride ((Ph₃P)AuCl), 98%

Part B

CAS : 14243-64-2

Molecular Formula : [(C₆H₅)₃P]AuCl

Molecular Weight : 494.71

General Information

Storage 8 to 25°C (Cool & Dry Area)

HSN Code

1 Gms

28439019 (GST 18%)

500 Mg

28439019 (GST 18%)

Available Packages

500 Mg

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com