



43384

## Hyaluronic Acid Sodium Salt - Type 5 (8-15 kD) extrapure, 90%

Part B

CAS : 9067-32-7

### Specifications

Appearance (Colour)	White
Appearance (Form)	Powder
Assay	min. 90%
Heavy metals	max. 0.001%

### Other Information

Application

Hyaluronic acid (HA) is an anionic linear glycosaminoglycan composed of monomers of glucuronic acid and N-acetylglucosamine. HA is used in many applications including, cell culture, tissue engineering, and the study of wound healing.

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Mg	39139090 (GST 18%)

### Available Packages

100 Mg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com