



43786

Maltotriose (DP3) extrapure, 90%

Part B

CAS : 1109-28-0

Molecular Formula : C₁₈H₃₂O₁₆

Molecular Weight : 504.44

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay	min. 90%
Specific rotation [?]	+164° to +166° (C=3 , Water)
Water (KF)	max. 5%

Other Information

Application

Maltotriose is a (1-4)-linked glucose oligomer that serves as the inducer of the maltose regulon of Escherichia coli.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	29400000 (GST 18%)
250 Mg	29400000 (GST 18%)
5 Gms	29400000 (GST 18%)

Available Packages

250 Mg

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com