



43800

Silver Hexafluoroantimonate(V) extrapure, 98%

Part B

CAS : 26042-64-8

Molecular Formula : AgSbF6

Molecular Weight : 343.62

Specifications

Appearance (Colour)	White to off-white to pale grey
Appearance (Form)	Crystals or powder
Assay (T)	min. 98%

Other Information

Application

Generates stable cation radical salts, is used as an acidic catalyst and to study the effect of ligands on network structure toward the goal of engineering novel materials.

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	28432900 (GST 18%)
5 Gms	28432900 (GST 18%)

Available Packages

1 Gms

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com