



43819

## Sodium Cobaltinitrite extrapure AR, ACS

Part A

CAS : 13600-98-1

Molecular Formula : Na<sub>3</sub>Co(NO<sub>2</sub>)<sub>6</sub>

Molecular Weight : 403.94

### Specifications

Appearance (Colour)	Orange-yellow to Brownish
Appearance (Form)	Crystalline powder
Insoluble matter	max. 0.02%
Chloride (Cl)	max. 0.005%
Sulphate (SO <sub>4</sub> )	max. 0.01%
Iron (Fe)	max. 0.002%
Suitability for determination of potassium	Passes test

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1479 , IMCO Class No.:5.1 , Packing Group:II
HSN Code	
100 Gms	28429090 (GST 18%)
25 Gms	28429090 (GST 18%)
500 Gms	28429090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

### Available Packages

25 Gms  
100 Gms  
500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com