



44000

Kohn Two Tube Medium No. 1 Base

Part D

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 8 to 25°C (Cool & Dry Area) |
| Shelf Life | 36 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 100 Gms | 38210000 (GST 18%) |
| 500 Gms | 38210000 (GST 18%) |

Available Packages

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com