



44111

Menadione Sodium Bisulphite Trihydrate W/S
(Vitamin K3 Sodium Bisulphite Trihydrate) pure,
95%

Part B

CAS : 6147-37-1

Molecular Formula : C₁₁H₈O₂.NaHSO₃.3H₂O

Molecular Weight : 330.26

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Assay (T)	min. 95%

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29362990 (GST 18%)
25 Gms	29362990 (GST 18%)
500 Gms	29362990 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com