



44170

Fullerene C70, 99%

Part C

Specifications

Appearance (Colour)	Black
Appearance (Form)	powder
Assay	min. 99%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	29420090 (GST 18%)
100 Mg	29420090 (GST 18%)
500 Mg	29420090 (GST 18%)

Available Packages

100 Mg
500 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com