



44369

N,N-Diethylaniline extrapure AR, 99%

Part A

CAS : 91-66-7

Molecular Formula : C₁₀H₁₅N

Molecular Weight : 149.24

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.933-0.935
Refractive Index (20°C)	1.541-1.542
Non Volatile Matter	max. 0.01%
Chloride (Cl)	max. 0.0005%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.001%
Ethylaniline+aniline	max. 1%
Suitability as a reagent for Zinc (Zn)	Passes test
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2432 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
2500 ml	29214190 (GST 18%)
500 ml	29214190 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com