



44445

5-Sulphosalicylic Acid extrapure, 99%

Part A

CAS : 5965-83-3

Molecular Formula : C7H6O6S.2H2O

Molecular Weight : 254.21

Specifications

Appearance (Colour)	White to off-white to grey
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Assay	min. 99%
Chloride (Cl)	max. 0.01%
Iron (Fe)	max. 0.001%
Heavy Metals (Pb)	max. 0.01%
Water (KF)	14 - 15%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2585 , IMCO Class No.:8 , Packing Group:III
HSN Code	
250 Gms	29182190 (GST 18%)
500 Gms	29182190 (GST 18%)
Type of Packing	
250 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

250 Gms

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com