



44683

Ammonium Acetate extrapure AR, 98%

Part A

CAS : 631-61-8

Molecular Formula : C2H7NO2

Molecular Weight : 77.08

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline hygroscopic powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (HCl titration)	min. 98%
Insoluble matter	max. 0.005%
pH (5% aq. solution)	6.5 - 7.5
Residue after ignition	max. 0.01%
Chloride (Cl)	max. 0.0005%
Sulphate (SO4)	max. 0.001%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.0005%
Nitrate (NO3)	max. 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
500 Gms	29152990 (GST 18%)

Available Packages

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.
608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com