



44826

L-Theanine extrapure, 99%

Part B

CAS : 3081-61-6

Molecular Formula : C₇H₁₅N₂O₃

Molecular Weight : 175.20

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Assay	min 99%
Specific rotation [?]	+7.7° to +8.5° (C=1, H ₂ O)
Loss on drying	max 0.2%
Chloride (Cl)	max 0.002%
Sulphate (SO ₄)	max 0.001%
Heavy Metals (Pb)	max 0.001%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29224290 (GST 18%)
5 Gms	29224290 (GST 18%)

Available Packages

5 Gms

25 Gms

100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com