



45057

## Pectin pure

Part A

CAS : 9000-69-5

Molecular Formula : (C<sub>7</sub>H<sub>20</sub>O<sub>6</sub>)<sub>130-260</sub>

Molecular Weight : 25000 to 50000

### Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Amorphous powder
Identification Test	Passes test
Loss on drying	max. 10%
Methoxyl content	min. 7%
Galacturonic acid	min. 80%
Starch	Passes test

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kg	13022000 (GST 18%)
100 Gms	13022000 (GST 18%)
500 Gms	13022000 (GST 18%)
Type of Packing	
1 Kg	Plastic Vial
100 Gms	Plastic Bottle
500 Gms	Plastic Bottle

### Available Packages

100 Gms

500 Gms

1 Kg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com