



45273

## Isopropylboronic Acid extrapure, 97%

Part B

CAS : 80041-89-0

Molecular Formula : C<sub>3</sub>H<sub>9</sub>BO<sub>2</sub>

Molecular Weight : 88.00

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline compound
Assay	min. 97%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
1 Gms	29319090 (GST 18%)
25 Gms	29319090 (GST 18%)
5 Gms	29319090 (GST 18%)

### Available Packages

1 Gms

---

5 Gms

---

25 Gms

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com