



46178

1-Butyl-1-Methylpyrrolidinium
Bis(trifluoromethylsulfonyl)imide (BMP TFSI)
extrapure, 98%

Part B

CAS : 223437-11-4

Molecular Formula : C₁₁H₂₀F₆N₂O₄S₂

Molecular Weight : 422.41

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay	min. 98%
Density (g/ml) @ 20°C	1.3980 to 1.4030 g/ml
Refractive Index (20°C)	1.4230 to 1.4270

Other Information

Application	For Catalysis and Nanotechnology.
-------------	-----------------------------------

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Gms	29339990 (GST 18%)
5 Gms	29339990 (GST 18%)

Available Packages

5 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com