



46220

DL-Dithiothreitol (DTT) ExiPlus, Multi-Compendial, 98%

Part B

CAS : 3483-12-3

Molecular Formula : C4H10O2S2

Molecular Weight : 154.24

Specifications

ExiPlus	meets compendial specs of BP
Appearance (Colour)	White
Appearance (Form)	Crystalline hygroscopic powder
Solubility (Turbidity) 5% aq. solution	Clear
Solubility (Colour) 5% aq. solution	Colourless
Assay (Iodometric)	min. 98%
Melting Point	41 - 43°C

Other Information

Application

Suitable for pharmaceutical production applications in bulk processing. For a reducing agent superior to DTT see Tris(2-carboxyethyl) phosphine hydrochloride(TCEP).

General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 Gms	29309099 (GST 18%)
Type of Packing	
5 Gms	Glass Bottle

Available Packages

5 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com